

# MIM Storage Server

## Comparison Chart

Included Free with  
MIM® Software

MIM Storage Server

### Features Included Free with MIM® Software

Fully Integrated Multi-modality Storage	X	X
No Additional Licensing for Adding Modalities	X	X
No Additional Licensing for 3rd Party Workstations	X	X
SSL Secure Network Transmission - 1024 bit	X	X
Remote MIM & MIMviewer™ Patient List Push/Pull	X	X
Advanced DICOM Reporting	X	X
DICOM Query/Retrieve Server	X	X
DICOM Store Server	X	X
DICOM Storage Commit	X	X
DICOM Printing	X	X
Age-based Auto-delete of Patient Studies	X	X

### MIM Storage Server

Role-based Security (Read/Write Privileges)	X
Integration with LDAP/Active Directory	X
Central MIM User & User Groups	X
Secure Networked User Login	X
Faster Query Performance Independent of Database Size	X
Automatic Patient Study Routing	X
Automatic Patient Study Retrieval from 3rd Party PACS	X
Per-user Configuration - Migrates to Remote Workstation	X
Age-based Auto-archive of Patient Studies	X
Web-based System Administration	Q4 2008
Configurable Intelligent Auto-Delete	Q4 2008
Attach Encapsulated User Documents in DICOM Format	Q4 2008
Archive Encapsulated User Documents in DICOM Format	Q4 2008